


Features

- Construction**— Carbon steel with non-skid surface protects animals and handlers from accident slips and falls. Adjustable foot pads allow scale to work effectively on uneven floors.
- Overload Capacity**— Capable of withstanding 1,500 pounds of weight without damaging the scale.
- Easy Assembly**— No installation costs. No moving parts that require maintenance or adjustment.
- Stainless Steel Junction Box**— allows easy connection to one of three Salter Brecknell indicators.
- Weight Sensor**— The Vet Deck uses a patented loadcell technology with proven performance in the animal industry for more than 30 years.




Specifications

- Capacity**— 1,000 lbs
- Accuracy**— 0.5%
- Dimensions**—
Scale: 50" (L) x 20" (W) x 1.75" (H)
Shipping Carton: 52" (L) x 21" (W) x 4" (H)
- Scale Weight**— 55 pounds
- Shipping Weight**— 62 pounds



- Indicator Option**— The Model S200E provides 5,000 divisions of certified accuracy (1000 x 0.2 lb.) with a bright LED display.



ONE DIGITAL SCALES.com

Your First Choice
for Industrial, Home and Specialty Scales

800-620-7811

www.onedigitalscales.com